




# LCMS- Report

LCMS  
Acquired by : SD  
Sample Name : M11151-IR808-03  
Method File : 3.lcm

Method info. : C18 4.6\*50mm 5um 2.000ml/min 2.6min Column Temperature: 40 °C  
Gradient: 5% B hold for 0.2min, increase to 95 % B within 1.40 min, hold at 95 % B for 0.9 min, then back to 5% B within 0.01 min  
Pump A : 0.03% TFA in H2O  
Pump B : 0.03% TFA in ACN

## Chromatogram



## Integration Results for DAD

PDA Ch1 254nm


| Peak# | Ret. Time | Height | Area   | Area%   |
|-------|-----------|--------|--------|---------|
| 1     | 1.948     | 88961  | 229895 | 100.000 |
| Total |           | 88961  | 229895 | 100.000 |

PDA Ch2 214nm

| Peak# | Ret. Time | Height | Area   | Area%   |
|-------|-----------|--------|--------|---------|
| 1     | 1.948     | 291126 | 759845 | 100.000 |
| Total |           | 291126 | 759845 | 100.000 |

### Mass Spectrum

Rt=1.947min Positive



Rt=2.058min Positive

