


<Chromatogram>


M11254-TUG-1375-01



Peak Table

PDA Ch1 254nm

Peak#	Ret. Time	Height	Height%	Area	Area%
1	0.503	1089	0.027	725	0.020
2	0.730	45390	1.123	26646	0.744
3	0.753	3995112	98.850	3552377	99.235
Total		4041590	100.000	3579748	100.000




Peak Table

AD2

Peak#	Ret. Time	Height	Height%	Area	Area%
1	0.738	1407578	100.000	4903853	100.000
Total		1407578	100.000	4903853	100.000

Segment#1 (x1,000,000)
Intensity

MS Chromatogram



Line#:1 R.Time:0.742(Scan#:61)

MassPeaks:814

Spectrum Mode:Single 0.742(61) BasePeak:443.1

BG Mode:None Segment 1 - Event 1

Mass Spectrum

