



Agilent Technologies

Sample name: M11262-MHP-01

Instrument: A010

Inj. volume: 4

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.757		3.8	1.2	0.049
6.745		0.8	0.3	0.011
7.003		3.0	1.0	0.038
7.610		7716.8	2173.4	99.775
8.223		1.6	0.5	0.021
12.316		8.2	2.1	0.107