




Sample name: M11346-ZN-c3-01
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
6.153		0.45	5283.7	1451.6	99.62
6.560		0.31	17.4	2.8	0.33
6.805		0.08	0.1	0.1	0.00
6.964		0.08	0.2	0.1	0.00
7.181		0.17	1.8	0.4	0.03
8.081		0.16	0.7	0.1	0.01