



Sample name: M11368-STM2457-01

Injection date:

Instrument: A012

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
4.111		20.7	4.2	0.151
4.245		19.7	6.7	0.143
4.494		5.4	1.1	0.039
4.778		16.4	3.8	0.119
5.182		13560.1	2761.7	98.446
5.374		36.3	6.7	0.263
5.821		100.5	27.0	0.729
6.094		6.4	1.9	0.046
7.513		8.6	2.5	0.062