



Agilent Technologies

Sample name: M11540-Propranolol-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.352	MB m	0.13	5.21	1.57	0.04
6.681	MM m	0.08	1.29	0.52	0.01
6.772	BM m	0.06	1.19	0.66	0.01
6.863	MM m	0.45	13817.97	3054.20	99.94