


Sample Name : M11553-WS3-01  
Acq. Operator : 9  
Spec. Reported : Integration

Location : P2-A-06  
Inj : 0  
Inj Volume : 5 ul



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
3.04	0.06	14.40	3.60	0.18
3.19	0.05	28.51	9.20	0.36
3.63	0.07	32.51	6.13	0.41
4.08	0.07	7809.41	1677.68	98.69
5.00	0.05	28.47	8.95	0.36