



Agilent Technologies

Sample name: M13561-C25-140-02

Instrument: A012

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
7.304		0.23	2.7	0.7	0.03
8.294		0.65	9985.1	2626.8	99.92
9.073		0.19	5.1	1.2	0.05