



Sample name: M13996-Pitolisant hydrochloride-02

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak	Area	Height	Area%
6.304		1.7	0.6	0.009
6.485		8.4	2.9	0.044
6.998		1.2	0.3	0.006
7.075		2.6	1.0	0.014
7.205		19.3	5.8	0.101
7.340		1.2	0.4	0.006
7.516		0.7	0.2	0.004
7.675		18964.4	2486.0	99.716
8.247		10.3	1.9	0.054
9.412		8.8	1.4	0.046