




Agilent Technologies

Sample name: M14034-Almorexant hydrochloride-02 Injection date:
Instrument: A010 Location:
Inj. volume: 2 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
4.945	1.43	0.54	0.01
6.575	4.64	2.00	0.04
6.634	1.77	0.80	0.02
6.909	11190.01	2239.95	99.80
7.344	8.61	1.53	0.08
9.194	5.59	0.90	0.05