


=====
Acq.Operator : System
Acq.Instrument : LC1100

Location : 68
Inj : 1
Actual Inj Volume : 3.000 µl



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier : 1.0000
Dilution : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=210,4 Ref=off

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	6.579	BV R	0.0560	7922.83643	2200.96826	99.2319
2	6.932	VB E	0.0669	15.11752	3.33569	0.1893
3	7.153	BV	0.0521	11.90998	3.46741	0.1492
4	7.234	VB	0.0572	34.30262	9.27050	0.4296

Totals : 7984.16653 2217.04186

=====
*** End of Report ***