




Sample name: M14158-JNK-IN-7-01
Instrument: A012
Inj. volume: 3
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
4.628	18.05	5.38	0.23
4.741	4.54	1.54	0.06
4.853	5.10	2.20	0.07
4.912	1.31	0.83	0.02
5.014	7559.02	1809.14	98.21
5.346	10.29	2.47	0.13
5.621	6.01	1.69	0.08
5.789	12.59	2.82	0.16
6.671	8.57	2.68	0.11
6.802	71.30	16.24	0.93