




Sample name: M14435-PF429242 dihydrochloride-01 **Injection date:**
Instrument: A012 **Location:**
Inj. volume: 3 **Type:**
Acq. method: HPLC **Sample amount:**
Processing method:
Manually modified:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.135		23.6	7.8	0.436
5.254		5338.9	1555.1	98.823
5.744		7.1	1.5	0.132
5.857		22.7	7.3	0.420
5.957		2.5	0.9	0.046
6.117		0.7	0.2	0.012
6.212		5.3	1.4	0.098
6.283		1.8	0.5	0.033