



Agilent Technologies

Sample name: M14523-DASA-58-01

Instrument: A010

Inj. volume: 2.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
6.176	3.66	1.26	0.04
6.814	9.01	2.57	0.10
7.355	8593.83	2206.39	99.75
7.853	2.06	0.64	0.02
8.000	1.40	0.62	0.02
8.127	5.48	1.91	0.06