



Injection date:

Sample name: M14629-Elacestrant dihydrochloride-02

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.962		2.1	0.5	0.028
6.060		3.2	1.3	0.042
6.196		7425.8	2207.5	99.508
6.431		3.2	1.0	0.043
6.572		1.8	0.7	0.023
6.698		18.9	6.4	0.253
6.853		7.5	2.6	0.101