



Agilent Technologies

Sample name: M14645-4-Azido-L-phenylalanine-01

Injection date:

Instrument: A010

Location:

Inj. volume: 5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak	Width [min]	Area	Height	Area%
4.651		0.15	3.4	1.0	0.04
4.806		0.14	3.0	0.9	0.04
5.287		0.58	8188.0	2545.7	98.13
6.808		0.21	47.2	14.7	0.57
7.394		0.28	12.1	2.5	0.14
8.644		0.32	90.7	20.6	1.09