




Sample Name : M15089-PBIT-01
Acq. Operator : 11
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%-95% B in 1.5min,95%B for 1.5 min
 Flow Rate:2ml/min
 Column:Sunfire ,4.6*50mm,3.5um
 Oven Temperature: 50 C
 MASS Range:120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
2.16	0.02	3884.55	2812.65	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.16	0.02	3024.12	2945.47	100.00

Ret. Time: 2.17

<<<< POSITIVE SPECTRA >>>>

