



Sample name: M15105-SR15006-02

Instrument: A012

Inj. volume: 10

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1F,Sig=260,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.625	MM m	0.21	4.49	0.89	0.03
6.882	BB	0.14	37.82	12.14	0.23
7.001	BB	0.10	8.90	3.30	0.05
7.158	BB	0.11	15.04	5.04	0.09
7.325	BB	0.44	16342.52	3367.45	99.16
7.704	BB	0.10	3.66	1.35	0.02
7.959	BB	0.25	51.56	15.62	0.31
8.261	BB	0.23	7.10	1.22	0.04
8.465	MM m	0.14	1.79	0.49	0.01
8.724	MM m	0.18	1.62	0.35	0.01
8.959	MM m	0.24	6.35	1.08	0.04