



Agilent Technologies

Sample name: M1663-AZD1152-HQPA-01

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.140	BB	0.16	5.25	1.54	0.06
5.346	BB	0.12	8.46	2.96	0.09
5.438	BB	0.08	5.10	2.07	0.06
5.598	BV	0.53	8884.06	2680.03	99.09
6.098	BB	0.15	18.47	4.26	0.21
6.254	BB	0.09	4.69	1.97	0.05
6.337	BB	0.10	5.95	2.23	0.07
6.866	BB	0.19	27.86	7.91	0.31
7.528	BV	0.14	5.48	1.52	0.06