



Sample name: M1673-Nintedanib (BIBF1120)-02

Injection date:

Instrument: A010

Location:

Inj. volume: 1


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
6.636	1.15	1.62	0.02
6.703	5582.36	1704.46	99.42
6.816	24.27	6.49	0.43
8.239	6.92	1.68	0.12