



Sample name: M1675-BIRB796-02

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.425	BB	0.14	63.81	21.59	0.53
5.544	MM m	0.11	2.27	0.75	0.02
6.356	BV	0.13	4.82	1.28	0.04
6.638	MM m	0.12	2.04	0.64	0.02
6.805	MM m	0.14	2.05	0.54	0.02
7.004	MM m	0.10	2.08	0.69	0.02
7.118	BB	0.13	6.62	1.50	0.06
7.315	BB	0.20	21.14	5.49	0.18
7.567	BV	0.46	11884.75	2850.30	99.00
9.183	BB	0.22	7.14	1.19	0.06
9.876	MB m	0.19	7.71	1.77	0.06