



Current Data Parameters
Name DCC-2036-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date 20210531
Time 9.27
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30 TD
65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 362
DW 60.400 usec
DE 8.00 usec
TE 297.2 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SF01 400.1324710 MHz

F2 - Processing parameters
SI 32768
SF 400.1300031 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

