


Sample Name : Entecavir Monohydrate-01

Method Info : Mobile Phase: A: water(10mMNH₄HCO₃) B: ACN
Gradient: 5%-95% B in 1.5min
Flow Rate: 1.8ml/min
Column:XBridge C18,4.6x50mm,3.5um
Oven Temperature: 50 C
UV214,MASS103-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.70	0.04	1384.83	525.51	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
0.70	0.04	1206.63	458.62	100.00

Ret. Time: 0.71 <<<< POSITIVE SPECTRA >>>>

