




Sample name: M1746-DAPT-08  
Instrument: A010  
Inj. volume: 2  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
8.360	8.53	2.52	0.06
8.803	34.48	6.61	0.23
10.005	14923.81	2150.96	99.50
10.157	5.65	1.48	0.04
10.644	26.74	5.29	0.18