



Agilent Technologies

Sample name: M1768-Rapamycin-10

Instrument: A010

Inj. volume: 1

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1G,Sig=280.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.331	BB	0.55	9.48	1.26	0.10
9.728	MM m	0.23	0.69	0.13	0.01
10.059	MM m	0.19	0.82	0.15	0.01
10.281	MM m	0.30	2.74	0.32	0.03
10.561	MM m	0.26	6.57	1.05	0.07
10.779	MM m	0.23	1.17	0.15	0.01
12.193	MB m	0.17	5.02	1.11	0.05
12.360	BM m	0.21	12.78	2.23	0.13
12.622	MM m	0.25	4.09	0.53	0.04
13.186	BB	1.84	9497.14	532.46	99.55