




Sample name: M1771-PTC124-01
Instrument: 3
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
6.471	2.01	0.80	0.02
9.052	1.26	0.36	0.01
9.389	9720.00	2032.78	99.41
9.751	6.32	1.78	0.06
9.905	0.96	0.37	0.01
10.005	3.76	1.18	0.04
12.161	43.66	9.78	0.45
Sum	9777.96		