




Sample name: M1773-PI-103-03
Instrument: A010
Inj. volume: 4
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A, Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
6.800		5910.9	1786.3	98.382
7.311		11.3	3.3	0.189
8.253		4.5	1.5	0.074
8.455		49.5	10.2	0.824
8.525		13.1	4.3	0.217
11.610		5.5	1.4	0.091
11.911		4.2	1.1	0.069
12.029		4.1	0.9	0.068
12.433		5.1	1.0	0.085