



Agilent Technologies

Sample name: M1774-PLX-4720-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
6.623	8.06	2.74	0.11
7.862	3.63	1.29	0.05
8.012	4.59	1.53	0.06
8.631	27.12	8.22	0.38
8.846	7022.12	1805.89	98.77
9.027	36.76	8.89	0.52
10.127	7.60	2.10	0.11