

Sample Name : M1776-Rivaroxaban-01-DMSO

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5% increase to 95%B within 1.3 min, 95%B for 1.7 min.


Flow Rate: 2.0ml/min

Column:Sunfire C18,4.6*50mm,3.5um


Column Temperature: 50 C

Detection:UV(214,4nm) and MS(ESI,Pos mode,50 to 700 amu)


DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.57	0.02	526.14	450.05	97.82

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.57	0.02	992.60	887.27	100.00

Ret. Time: 1.59

<<<< POSITIVE SPECTRA >>>>

