



Agilent Technologies

Sample name: M1801-KU-55933-01

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.265	MM m	0.36	8.24	0.92	0.09
6.717	MM m	0.09	1.12	0.35	0.01
6.925	MM m	0.09	1.41	0.46	0.02
7.116	MB m	0.17	9.84	1.93	0.10
7.285	BB	0.16	8.09	2.09	0.09
8.448	BM m	0.21	9287.09	2524.08	98.61
8.575	MM m	0.11	92.43	28.00	0.98
9.885	MM m	0.13	1.82	0.44	0.02
11.544	MM m	0.16	1.72	0.35	0.02
12.867	MM m	0.15	6.26	1.45	0.07