


File D:\DATA\
 Sample Name : M1811-ARQ-197-01
 Acq. Method : D: HPLC -15MIN.M
 Analysis Method : D: HPLC -15MIN.M
 Sample Info : Easy-Access Method:
 Method Info : Column:Xbridge C18 (4.6x150 mm, 3.5um)
 Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)
 Elution program: Gradient from 5 to 95% of B in 9min at 1ml/min
 Temperature: 50°C
 Detection: UV (214 , 4 nm) and MS (ESI, Neg mode , 70 to 900 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
8.29	0.04	5349.18	2435.50	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
8.29	0.03	626.32	298.58	100.00