




Agilent Technologies

Sample name: M1817-Y-27632 dihydrochloride-30-1 Injection date:
Instrument: A010 Location:
Inj. volume: 2 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1E, Sig=260.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
1.526	BM m	0.66	10476.42	2799.91	99.96
1.975	MM m	0.29	3.67	0.42	0.04