



Agilent Technologies

Sample name: M1825-PCI-32765-03

Instrument: A012

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.119	MB m	0.13	8.87	2.75	0.11
7.509	BB	0.34	8218.03	2324.31	99.83
7.986	MB m	0.13	4.76	1.34	0.06