

M1828-Tandutinib-01-LCMS

Sample Information

LCMS


Acquired by : System Administrator
Sample Name : M1828-Tandutinib-01
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/7/17 10:12:55
Injection Volume : 3
Vial : 25


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.731	21512	415580	100.000
总计		21512	415580	100.000

PDA Ch2 214nm


Peak#	Ret. Time	Height	Area	Area%
总计				

Mass Spectrum


Rt=0.394min 正




Rt=0.557min 负




Rt=0.912min 正




Rt=1.165min 负




Rt=1.384min 正




Rt=1.637min 负




Rt=1.699min 正



Rt=2.110min 负



Rt=2.267min 負



Rt=2.352min 正

