


Sample Name : M1837-MK-2206 2HCl-06
Acq. Operator : 9
Spec. Reported : Integration

Location : P1-C-05
Inj : 0
Inj Volume : 5 ul



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
3.04	0.08	17.95	3.66	0.18
3.14	0.06	29.42	6.50	0.29
3.37	0.05	34.85	10.59	0.35
3.48	0.06	9931.21	2489.10	99.18