



Agilent Technologies

Sample name: M1850-WZ4002-09

Instrument: A010

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Peak	Area%
6.754	13.66	3.98		0.18
7.338	7514.06	2113.05		99.35
7.579	3.22	1.39		0.04
7.649	9.93	3.78		0.13
7.906	7.03	2.27		0.09
8.485	15.61	3.86		0.21