




Agilent Technologies

**Sample name:** M1861-Silmitasertib (CX-4945)-04-2 **Injection date:**  
**Instrument:** A010 **Location:**  
**Inj. volume:** 8 **Type:**  
**Acq. method:** HPLC **Sample amount:**  
**Processing method:**  
**Manually modified:**



**Signal:** DAD1E,Sig=260,4 Ref=360,100

| RT [min] | Type | Width [min] | Area    | Height  | Area% |
|----------|------|-------------|---------|---------|-------|
| 7.300    | BM m | 0.99        | 7230.61 | 1205.96 | 99.78 |
| 10.687   | MB m | 0.16        | 15.82   | 3.70    | 0.22  |