



Agilent Technologies

Sample name: M1862-GDC-0068-03

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Peak	Area%
5.118	3.80	1.45		0.06
5.706	22.91	7.90		0.33
6.704	2.38	0.98		0.03
8.742	6856.81	1153.40		99.34
9.125	16.43	4.52		0.24