



Agilent Technologies

Sample name: M1875-PD173074-02

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
5.765	MM m	0.10	86.29	31.04	1.28
6.827	BM m	0.29	6635.16	1613.69	98.16
7.066	MM m	0.08	12.80	3.88	0.19
7.581	BM m	0.09	6.83	2.65	0.10
8.076	MB m	0.10	6.02	2.05	0.09
9.265	MM m	0.12	12.38	3.64	0.18