



Sample name: M1924-Linezolid-06

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
4.468		0.10	0.3	0.1	0.00
5.865		0.51	8337.4	2565.5	99.77
6.248		0.13	2.4	0.5	0.03
6.408		0.11	0.6	0.1	0.01
7.390		0.11	2.1	0.7	0.02
8.211		0.26	13.6	3.3	0.16