



Sample name: M1925-LY294002-08

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.976	MM m	0.07	0.66	0.24	0.00
7.740	MM m	0.36	13973.95	2293.42	99.95
8.180	BB	0.14	3.77	1.08	0.03
8.595	MM m	0.11	2.16	0.61	0.02
8.733	MM m	0.10	0.69	0.19	0.00