

LCMS- Report




Sample Information

LCMS
Acquired by : SS
Sample Name : M1940-Docetaxel-03
Method File : B.lcm

Method info. : C18 4.6*50mm 5um 2.000ml/min 2.6min Column Temperature: 40 °C
Gradient: 5% B hold for 0.2min, increase to 95 % B within 1.40 min, hold at 95 % B for 0.9 min, then back to 5% B within 0.01 min

Pump A : 0.03% TFA in H2O
Pump B : 0.03% TFA in ACN

Chromatogram



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.571	14049	37164	1.264
2	1.744	1068286	2902277	98.736
Total		1082336	2939441	100.000

PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.746	2000070	9710254	100.000
Total		2000070	9710254	100.000

Mass Spectrum

Rt=1.492min Positive



Rt=1.743min Positive

