



Sample name: M20182-WAY-600704-01

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.333	MM m	0.08	0.36	0.12	0.00
6.538	MM m	0.14	2.95	0.76	0.02
7.111	MM m	0.41	16301.45	3363.95	99.94
7.662	MM m	0.09	0.53	0.18	0.00
8.867	MM m	0.17	2.44	0.59	0.01
9.653	MM m	0.12	0.68	0.18	0.00
9.784	MM m	0.13	2.65	0.70	0.02