



Agilent Technologies

Sample name: M20602-UM729-02

Instrument: A010

Inj. volume: 2.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1F,Sig=260,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.587	BM m	0.08	33.51	16.84	0.31
6.670	MM m	0.23	10737.50	2977.67	99.16
6.851	MM m	0.14	53.08	13.20	0.49
7.285	BM m	0.14	4.33	1.04	0.04