



Agilent Technologies

Sample name: M20693-S-Gboxin-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
6.421	10.09	4.22	0.14
6.984	6.34	2.55	0.09
7.188	9.70	3.78	0.13
11.132	7363.29	786.86	99.40
12.054	18.04	5.02	0.24