



Sample name: M2081-LY2109761-07

Injection date:

Instrument: 1.5

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
4.805		0.05	3919.8	1287.3	98.19
5.089		0.04	5.2	2.1	0.13
5.497		0.04	61.5	25.7	1.54
5.661		0.04	5.8	2.4	0.14
		Sum	3992.18		