



Sample name: M2087-HC-030031-02

Injection date:

Instrument: 2

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
7.231		0.20	5.5	1.3	0.04
7.766		0.42	12495.4	3173.7	99.94
8.385		0.06	1.6	0.5	0.01
		Sum	12502.52		