


Sample name: M20874-Elamipretide-01
 Description: Mobile phase: Buffer A: 0.1% TFA in H₂O
 Buffer B: 0.1% TFA in ACN
 Gradient: 15-25% B in 20 min, Flow rate: 1 ml/min
 Column: Unitary C18, 4.6*150 mm, 5 µm 100 Å
 Instrument: Agilent HPLC 1100

Instrument: 1100

Location: 91
 Injection: 1 of 1
 Injection volume: 5.000
 Acq. operator: SYSTEM



Signal: DAD1 A, Sig=220, 4 Ref=360, 100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
9.350	VV R	0.3213	12649.3359	556.3600	99.9619	
13.263	MM	0.1978	4.8256	0.4066	0.0381	
Sum			12654.1615			