

MASS SPECTROMETRY REPORT



Sample Description

Analyst: SHAO
Sample: M20874-Elamipretide-01
M.W.: 639.8
Lot. No.:

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters 2695

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H2O/50%ACN