



Sample name: M2089-Irinotecan Hydrochloride Trihydrate-01

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.640		1.5	0.6	0.012
5.753		16.4	4.8	0.128
5.902		10.7	3.0	0.083
6.032		12731.2	2615.8	99.445
6.339		8.4	2.9	0.066
6.441		15.7	5.4	0.123
6.537		8.3	2.8	0.065
7.118		0.9	0.3	0.007
7.292		1.9	0.6	0.015
7.431		2.7	0.8	0.021
7.714		4.5	1.3	0.035